

Setup Opal-RAD Workstation for viewing - Internet Explorer 10/9

NOTE: Internet Explorer must be used, all other web browsers are not compatible for full viewing*

1. Open Internet Explorer Browser
2. Type the location of the Server in the Address bar:
(ie. http://OPAL_QC/opalweb/login.aspx or <http://IP Address/opalweb/login.aspx>) and hit Enter



If you receive a prompt, select 'Turn on Intranet Settings'


You may skip step 3 if Intranet settings are enabled (enables for Server connection only)

Intranet settings are turned off by default.


Don't show this message again

Turn on Intranet settings

3. Enable 'Compatibility View'
 1. Look for the Compatibility View button . It will appear in right side of the Address bar
(If you don't see the button, there's no need to turn on Compatibility View)
 2. Click the Compatibility View button  to enable Compatibility View for the Server connection
4. Restore down the browser window (so it is not full screen)
5. Select 20/20 logo/icon to the left of the address bar, click and drag it onto the Desktop to create a shortcut



6. Login
7. You will be prompted to install '**OpalViewer**'. This is a safe program. Click OK to allow install
**If OpalViewer prompt does not appear: Click upper right, 'Download Viewer - 2.4.x.x'*



8. After OpalViewer is installed, you may be asked to allow '**OpalWebKicker**'. This is necessary for viewing software to run correctly

This webpage wants to run the following add-on: 'OpalWebKicker' from 'Viztek LLC (unverified publisher)'. [What's the risk?](#)

Allow

Click 'Allow'. The session might log you off, this is normal. Simply log back in

9. Your viewer is now complete

*mOpal is available for other browsers as well as mobile and smart devices. See website or contact 20/20 Imaging for details